

## CA-1000 Capillary Columns

- Polyethylene glycol nitroterephthalic acid (also known as FFAP)
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness and efficiency

### Temperature Limits (0.25µm film)

Minimum	40°C
Maximum Isothermal	250°C
Maximum Programmed	260°C

### Similar Phases

AT-1000, BP-21, CP-Wax 58CB, DB-FFAP, HP-FFAP, Nukol, Stabilwax-DA, USP G25, USP G36, VB-FFAP, ZB-FFAP

### Common Applications

Alcoholic beverages and free fatty acids

Length	ID	Film	Catalog Number	Price
15	0.25	0.25	1503605	\$325.00
30	0.25	0.25	1506605	\$505.00
60	0.25	0.25	1509605	\$875.00
15	0.25	1.00	1503607	\$325.00
30	0.25	1.00	1506607	\$505.00
60	0.25	1.00	1509607	\$875.00
15	0.32	0.25	1503705	\$335.00
30	0.32	0.25	1506705	\$525.00
60	0.32	0.25	1509705	\$950.00
60	0.32	0.30	1509721	\$950.00
15	0.32	1.00	1503707	\$335.00
30	0.32	1.00	1506707	\$525.00
60	0.32	1.00	1509707	\$950.00
15	0.53	0.50	1503907	\$350.00
30	0.53	0.50	1506907	\$600.00
60	0.53	0.50	1509907	\$975.00
15	0.53	1.00	1503907	\$350.00
30	0.53	1.00	1506907	\$600.00
60	0.53	1.00	1509907	\$975.00
30	0.53	1.20	1506919	\$600.00