

CA-1701 Capillary Columns

- 14% Cyanopropyl silicone
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness and efficiency

Temperature Limits (0.25 μ m film)

Minimum	-60°C
Maximum Isothermal	270°C
Maximum Programmed	280°C

Similar Phases

AT-1701, BP-10, CP Sil 19CB, DB-1701, OV-1701, RTx-1701, SPB-1701, USP G46, VB-1701, ZB-1701, 007-1701

Common Applications

Herbicides, pesticides, drugs and solvents

Length	ID	Film	Catalog Number	Price
15	0.25	0.25	0903605	\$285.00
30	0.25	0.25	0906605	\$485.00
60	0.25	0.25	0909605	\$850.00
15	0.25	1.00	0903607	\$285.00
30	0.25	1.00	0906607	\$485.00
60	0.25	1.00	0909607	\$850.00
15	0.32	0.25	0903705	\$305.00
30	0.32	0.25	0906705	\$505.00
60	0.32	0.25	0909705	\$910.00
15	0.32	1.00	0903707	\$305.00
30	0.32	1.00	0906707	\$505.00
60	0.32	1.00	0909707	\$910.00
15	0.53	1.00	0906907	\$325.00
30	0.53	1.00	0906907	\$575.00
60	0.53	1.00	0909907	\$950.00
30	0.53	1.20	0906919	\$575.00
15	0.53	1.50	0903908	\$325.00
30	0.53	1.50	0906908	\$575.00
60	0.53	1.50	0909908	\$950.00