

CA-1 Capillary Columns

- Methyl silicone
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness and efficiency

Temperature Limits (0.25 μ m film)

Minimum	-60°C
Maximum Isothermal	340°C
Maximum Programmed	350°C

Similar Phases

AT-1, BP-1, CP Sil 5CB, DB-1, Equity-1, HP-1, MDN-1, OV-1, RTx-1, SE-30, SPB-1, Ultra-1, USP G1, USP G2, USP G9, VB-1, VF-1, ZB-1, 007-1

Common Applications

Petroleum, petrochemical

Length	ID	Film	Catalog Number	Price
15	0.25	0.25	0103605	\$285.00
30	0.25	0.25	0106605	\$485.00
60	0.25	0.25	0109605	\$850.00
100	0.25	0.50	0111606	\$975.00
15	0.25	1.00	0103607	\$285.00
30	0.25	1.00	0106607	\$485.00
60	0.25	1.00	0109607	\$850.00
15	0.32	0.25	0103705	\$305.00
30	0.32	0.25	0106705	\$505.00
60	0.32	0.25	0109705	\$910.00
15	0.32	1.00	0103707	\$305.00
30	0.32	1.00	0106707	\$505.00
60	0.32	1.00	0109707	\$910.00
15	0.53	1.50	0103908	\$325.00
30	0.53	1.50	0106908	\$575.00
60	0.53	1.50	0109908	\$950.00
15	0.53	5.00	0103910	\$325.00
30	0.53	5.00	0106910	\$575.00
60	0.53	5.00	0109910	\$950.00
100	0.53	5.00	0111910	\$1,195.00