

CA-2560 Capillary Columns

- 100% biscyanopropyl polysiloxane
- Individually tested for inertness and efficiency

Temperature Limits (0.25 μ m film)

Minimum 20°C

Maximum Isothermal 250°C

Similar Phases

CP-sil 88, HP-88, RT-2560, SP-2560

Common Applications

Cis/trans FAMES

Length	ID	Film	Catalog Number	Price
30	0.25	0.20	3106604	\$460.00
60	0.25	0.20	3109604	\$800.00
100	0.25	0.20	3111604	\$1,275.00