

CA-502.2 Capillary Columns

- Proprietary phenyl methyl silicone
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness, efficiency and bleed

Temperature Limits (0.25µm film)

| | |
|--------------------|-------|
| Minimum | -20°C |
| Maximum Isothermal | 250°C |
| Maximum Programmed | 270°C |

Similar Phases

AT-502.2, DB-502.2, RTx-502.2

Common Applications

Purgeable volatiles

| Length | ID | Film | Catalog Number | Price |
|--------|------|------|----------------|------------|
| 30 | 0.25 | 1.40 | 0706613 | \$525.00 |
| 60 | 0.25 | 1.40 | 0709613 | \$865.00 |
| 30 | 0.32 | 1.80 | 0706714 | \$550.00 |
| 60 | 0.32 | 1.80 | 0709714 | \$935.00 |
| 30 | 0.53 | 3.00 | 0706909 | \$610.00 |
| 60 | 0.53 | 3.00 | 0709909 | \$975.00 |
| 105 | 0.53 | 3.00 | 0712909 | \$1,250.00 |