

CA-5 Capillary Columns

- 5% Phenyl methyl silicone
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness and efficiency

Temperature Limits (0.25µm film)

Minimum	-60°C
Maximum Isothermal	325°C
Maximum Programmed	340°C

Similar Phases

AT-5, BP-5, CP Sil 8CB, DB-5, Equity-5, HP-5, MDN-5, RTx-5, SE-54, SPB-5, Ultra-2, USP G27, USP G36, VB-5, VF-5, XTI-5, ZB-5, 007-5

Common Applications

Herbicides, pesticides, polyaromatic hydrocarbons, drugs and flavors

Length	ID	Film	Catalog Number	Price
15	0.25	0.25	0203605	\$285.00
30	0.25	0.25	0206605	\$485.00
60	0.25	0.25	0209605	\$850.00
15	0.25	1.00	0203607	\$285.00
30	0.25	1.00	0206607	\$485.00
60	0.25	1.00	0209607	\$850.00
15	0.32	0.25	0203705	\$305.00
30	0.32	0.25	0206705	\$505.00
60	0.32	0.25	0209705	\$910.00
15	0.32	1.00	0203707	\$305.00
30	0.32	1.00	0206707	\$505.00
60	0.32	1.00	0209707	\$910.00
15	0.53	1.50	0203908	\$325.00
30	0.53	1.50	0206908	\$575.00
60	0.53	1.50	0209908	\$950.00
15	0.53	5.00	0203910	\$325.00
30	0.53	5.00	0206910	\$575.00
60	0.53	5.00	0209910	\$950.00