

### CA-624 Capillary Columns

- Cyanopropyl silicone
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness and efficiency

### Temperature Limits (0.25µm film)

Minimum	-20°C
Maximum Isothermal	240°C
Maximum Programmed	260°C

### Similar Phases

AT-624, DB-624, RTx-624, ZB-624, G43

### Common Applications

EPA Method 624, purgeable organics from water

Length	ID	Film	Catalog Number	Price
30	32	1.80	0506714	\$575.00
60	0.32	1.80	0509714	\$950.00
30	0.53	3.00	0506909	\$610.00
60	0.53	3.00	0509909	\$950.00
75	0.53	3.00	0510909	\$1,175.00
100	0.53	3.00	0511909	\$1,250.00