

CA-Wax Capillary Columns

- Polyethylene glycol
- Fully crosslinked and surface bonded
- Solvent rinseable
- Individually tested for inertness and efficiency

Temperature Limits (0.25 μ m film)

Minimum	40°C
Maximum Isothermal	270°C
Maximum Programmed	280°C

Similar Phases

AT-Wax, BP-20, BTR-CW, CP-Wax 52CB, DB-Wax, DB-WAXetr, HP-INNOwax, Stabilwax, Supelcowax-10, USP G16, VB-Wax, 007-CW

Common Applications

Alcoholic beverages, foods, fragrances, flavors, fatty acid methyl esters and solvents

Length	ID	Film	Catalog Number	Price
15	0.25	0.25	1403605	\$325.00
30	0.25	0.25	1406605	\$505.00
60	0.25	0.25	1409605	\$875.00
15	0.25	1.00	1403607	\$325.00
30	0.25	1.00	1406607	\$505.00
60	0.25	1.00	1409607	\$875.00
15	0.32	0.25	1403705	\$335.00
30	0.32	0.25	1406705	\$525.00
60	0.32	0.25	1409705	\$950.00
15	0.32	1.00	1403707	\$335.00
30	0.32	1.00	1406707	\$525.00
60	0.32	1.00	1409707	\$950.00
15	0.53	0.50	1403907	\$350.00
30	0.53	0.50	1406907	\$600.00
60	0.53	0.50	1409907	\$975.00
15	0.53	1.00	1403907	\$350.00
30	0.53	1.00	1406907	\$600.00
60	0.53	1.00	1409907	\$975.00
30	0.53	1.20	1406919	\$600.00